

CUSTOMER SERVICE BULLETIN

Radio Shack® CUSTOMER SERVICE BULLETIN

400 Atrium—One Tandy Center, Fort Worth, Texas 76102

Bulletin No.	1568-β2
Description	Medical Office Systems
Model	Model III
Stock No.	26-1568 Version 1.β
Date	Dec. 31, 1981
Page 1 of	1

Problem

During edit transactions, if an account has two transactions and you tell it to edit line 3 (to add a line), after the insertion the line is found on the next account.

Solution

To correct this problem the following should be applied.

1. At TRSDOS READY, type BASIC.
2. Type LOAD"EDITOR".
3. Add the following line.
1185 IFCL>NLTHEN116β
4. Type SAVE"EDITOR" to save the change.
5. At TRSDOS READY, make a backup copy of the corrected diskette.